# LiNA Skin Hook™ (non-sterile)

Ref: SH-110-NS, SH-220-NS

## Intended use:

The LiNA Skin Hook™ is intended to retract the soft tissue from the operative surface.

## Patient population:

Patients intended to have surgical procedures which require optimal exposure.

## Product description:

The LiNA Skin Hook™ together with LiNA SeaStar™ is a complete, self-retaining system for retracting delicate soft tissue. By securing the LiNA Skin Hook™ on to the LiNA SeaStar™ the exposure of the incision is increased. The flexible and elastic hook can be adjusted during the procedure to create optimum tension.

## Precautions:

- For single use only. Do not reuse, reprocess or re-sterilize this device. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction.
- Use of the LiNA Skin Hook™ requires adequate training and experience in performing surgical procedures. Intended user profile is a trained surgeon or physician.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy. Be aware of the sharp edges.
- Please do not attach LiNA Skin Hook™ to home-made retractor rings.

## Warnings:

- Do not use past expiration date.
- Carefully inspect the device for any damage prior to use. Dispose of the device if bent, malformed or if breaks or cracks are found.
- Special attention is needed when installing and removing LiNA Skin Hook™ from the operative surface since device has sharp edges.

## Instruction for use:

1. Device is provided as a non-sterile therefore for devices packed in sterilization pouch before first use please use recommended sterilization process parameters:

- Product is designated for EtO sterilization.
- Product should not be cleaned before sterilization.

### Option A:

- Maximum load density 103kg/m<sup>3</sup>.
- Preconditioning: Temp. 35-45°C; Humidity 45-75%. Min. preconditioning time: 16 hours.
- Exposure phase:

Min. time for humidity: 00:30 min.; EtO dwell time: 04:59-05:10h;  
EtO concentration during gas dwell phase: 582-794 mg/l; temp. during dwell time: 40-50°C.

- Aeration:

Temp. 35-50°C; time of aeration: min. 24h.

### Option B:

- Maximum load density 161.7 kg/m<sup>3</sup>.
- Preconditioning: Temp. 40-50°C; Humidity 55-75%. Min. preconditioning time: 10 hours.
- Exposure phase:

Min. time for humidity: 15 min.; EtO dwell time: 04:50-05:00h;  
EtO concentration during gas dwell phase: 884,7-901,7 mg/l; temp. during dwell time: 40-50°C.

- Aeration:

Desorption via 20-30h in temp. 40-50°C and then quarantine stage via 6 days in room temperature.

## Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical ApS and the competent regulatory authority of the country in which the user and/or patient is established.

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