

LiNA Librata™

Innovation in Gynecology

LiNA Librata™

REF: LIB-1, LIB-3

INDICATION FOR USE:

The LiNA Librata™ is a thermal balloon ablation device intended to ablate the endometrial lining of the uterus in premenopausal women with menorrhagia (excessive uterine bleeding) due to benign causes for whom childbearing is complete.

CONTRAINDICATIONS:

The device is contraindicated for use in:

- A patient who is pregnant or who wants to become pregnant in the future. Pregnancies following ablation can be dangerous for both mother and fetus.
- A patient with any anatomic condition (e.g., history of previous classical cesarean sections or transmurall myomectomy) or pathologic condition (e.g., chronic immunosuppressive therapy) that could lead to weakening of the myometrium.
- A patient with known or suspected endometrial carcinoma (uterine cancer) or premalignant change of the endometrium, such as unresolved adenomatous hyperplasia.
- A patient with an intrauterine device (IUD) currently in place.
- A patient with active genital or urinary tract infection at the time of procedure (e.g., cervicitis, vaginitis, endometritis, salpingitis, or cystitis) or with active pelvic inflammatory disease (PID).
- Repeated ablation should not be performed using global ablation but only by hysteroscopic ablation by an experienced surgeon.

WARNINGS:

Failure to follow all instructions or to heed any warnings or precautions could result in serious patient injury.

- The LiNA Librata™ is provided STERILE. Carefully inspect the packaging for any damage prior to use. Do NOT attempt to use the device if sterile barrier is damaged. Do NOT use past expiration date. Do NOT use if the device is exposed to non-sterile surfaces before procedure.
- For single use only. Do NOT reuse, reprocess or re-sterilize the LiNA Librata™. Any reprocessing may impede the functions of this device.
- Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction and/or erroneous pathology specimen collection due to residual tissue in the LiNA Librata™.
- No modification of this equipment is allowed.
- The ambient temperature must not be higher than 30°C during the use of this equipment.
- The device is intended for use only in women who do not desire to bear children. There have been reports of women becoming pregnant following this procedure. Pregnancies after ablation can be dangerous for both mother and fetus.
- Endometrial ablation using the LiNA Librata™ is not a sterilization procedure. An effective form of contraception is required following LiNA Librata™ procedure.
- Endometrial ablation procedures using the LiNA Librata™ should be performed only by medical professionals who have experience in performing procedures within the uterine cavity, such as IUD insertion or dilation and curettage (D&C), and who have adequate training and familiarity with the LiNA Librata™.

- Do NOT perform same-day LiNA Librata™ procedures and hysteroscopic tubal occlusion/sterilization. Ablation may cause intrauterine synechiae, which can compromise (i.e., prevent) the 3-month confirmation test (HSG) for the tubal occlusion device. Women who have inadequate 3-month confirmation tests cannot rely on the tubal occlusion device for contraception.
- Patients who undergo endometrial ablation procedures who have previously undergone tubal ligation are at increased risk of developing post-ablation-tubal sterilization syndrome, which can require hysterectomy. This can occur as late as 10 years post-procedure.
- Uterine perforation can occur during any procedure in which the uterus is instrumented. Use caution not to perforate the uterine wall when sounding the uterus, dilating the cervix, curettaging the uterine lining or inserting the catheter. A hysteroscopy is recommended prior to LiNA Librata™ insertion to ensure that the uterus has not been perforated during sounding, dilation or curettage. If a uterine perforation is suspected after LiNA Librata™ insertion, THE PROCEDURE SHOULD BE TERMINATED IMMEDIATELY and hysteroscopy should be performed to assess the uterine cavity.
- If a perforation is present, and the procedure is not terminated, thermal injury to adjacent tissue may occur.
- Endometrial ablation procedures do not eliminate the potential for endometrial hyperplasia or adenocarcinoma of the endometrium and may mask the physician's ability to detect or make a diagnosis of such pathology.
- Once treatment is initiated in the uterus, patients are to be treated using 1 LiNA Librata™ device. LiNA Librata™ delivers an automated, 5 treatment cycle therapy. If the treatment is interrupted during any of the 5 cycles, DO NOT treat with any additional devices because of the potential for injury. Do not treat patients for more than one therapy cycle in a given treatment session because of the potential for transmural injury to the uterus or injury to adjacent viscera.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of LiNA Librata™. Otherwise, degradation of the performance of this equipment could result.
- The balloon must not come into contact with sharp edges.
- The distal end of the housing can reach temperatures of up to 48°C and should not be held.
- The balloon is very hot during and right after treatment.
- Do not remove the balloon catheter during the treatment process.
- Do not remove the battery during the treatment process.
- If a balloon leak is suspected during treatment, abort the procedure by following the instructions below in 'To Abort the Procedure'.
- If unexpected resistance is felt prior to the removal of the catheter pull the "EMERGENCY" strap on the back of the device until the word "STOP" appears. Wait 5 seconds and remove the catheter slowly from the uterus.
- Do not insert the balloon catheter into the patient if there is glycerin and/or air in the balloon.

LiNA Librata™

REF: LIB-1, LIB-3

POTENTIAL RISKS:

As with any procedure in which the uterus is instrumented certain risks exist. Hysteroscopy, and endometrial ablation are well-known procedures. Risks associated with the LiNA Librata™ device and/ or endometrial procedure, include but may not be limited to the following:

- Nausea
- Vomiting
- Urinary Tract Infection (UTI)
- Laceration
- Abdominal Pain
- Thermal injury to adjacent tissue, including cervix, vagina, vulva and/or perineum, fallopian tubes and adjacent organs (e.g. bowel burn)
- Vaginal infection
- Endometritis
- Hematometra
- Diarrhea
- Fever
- Headache
- Abdominal distension
- Perforation of uterus or creation of a false passage within the uterine wall
- Risks associated with hysteroscopy
- Infection or sepsis
- Complications leading to serious injury or death
- Myometritis / myonecrosis
- Cervical stenosis / intrauterine adhesions
- Pelvic inflammatory disease / pelvic abscess
- Vaginal discharge
- Abdominal cramping
- Vasovagal episodes
- Perioperative haemorrhage
- Transient urinary incontinence
- Voiding dysfunction

PRECAUTIONS:

- The safety and effectiveness of the LiNA Librata™ has not been fully evaluated in patients:
 - with uterine sounding >10 cm
 - with uterine sounding <6 cm
 - with bicornuate or septate uteri or previous endometrial resection/ablation
 - who are post-menopausal
- A patient with submucousal / intramural myomas such that the uterine cavity is significantly distorted.
- Uterine wall perforation or the creation of a false passage within the uterine wall can occur during any procedure in which the uterus is instrumented, especially in cases of severe anteverted, retroflexed, or a laterally displaced uterus.
- It is important that when the LiNA Librata™ balloon catheter is advanced in to the uterus that the indicator ring (set at the depth obtained during sounding) reaches the cervix without meeting resistance. If resistance is felt prior to the indicator ring reaching the cervix, the balloon catheter should be removed and a second sounding should be performed and the indicator ring adjusted accordingly.
- Use caution to ensure that the device is properly positioned in the uterine cavity, ensure that the catheter advances to the depth that was obtained during the sounding procedure.
- It has been reported that patients with a severe anteverted, retroflexed, or laterally displaced uterus are at an increased risk of uterine wall perforation

during any intrauterine manipulation. The clinician should use discretion in patient selection.

- History of repeat low transverse Cesarean section may increase the risk of patient injury due to thinning of the myometrium.
- Please note that portable and mobile RF communications equipment can affect the LiNA Librata™ device.
- Always have a backup device readily available for use.

TECHNICAL DESCRIPTION:

The LiNA Librata™ is a Type BF instrument, according to IEC 60601-1.

Technical data:

Essential performance	
Maximum glycerin chamber temperature	160 ±5 °C
Maximum treatment pressure	187 ±10 mmHg
Pressurized treatment time	126 seconds
Battery rating	
Nominal Voltage	14,8 V
Nominal capacity	2,9 Ah
Nominal Energy	42,92 Wh
Environmental conditions during transportation and storage	
Temperature	0-38°C
Air humidity	0-85% RH (non-condensing)
Atmospheric pressure	700-1060 hPa (0 to 3,030 m)
Environmental conditions during use	
Temperature	10-30°C
Air humidity	30-75% RH (non-condensing)
Atmospheric pressure	700-1060 hPa (0 to 3,030 m)

INSTRUCTIONS FOR USE:

Read the instruction for use prior to using this device!

A thinned endometrium may enhance the results of the LiNA Librata™ treatment. This can be accomplished by timing the menstrual cycle to the early proliferative phase, administering pre- treatment drugs such as danocrine, birth control pills, provera, etc. GnRH agonists should not be used to thin the endometrium. The optimum pre- treatment regimens have not been determined at this time.

To abort the procedure:

- Press the "ON" and "TREAT" buttons simultaneously ("ABORT")
- Observe that the "FINALIZED" light is flashing and that beeps are sounding.
- Wait 8 seconds to let the "ABORT" feature withdraw all the fluid from the balloon. The "FINALIZED" light will illuminate.

LiNA Librata™

REF: LIB-1, LIB-3

- D. Remove the catheter slowly from the uterus.
 - E. In the unlikely event that the "ABORT" function does not work, or if unexpected resistance is felt during removal of the catheter, pull the "EMERGENCY" strap on the back of the device until the word "STOP" appears. This will drain the fluid from the balloon back into device.
 - F. Wait 5 seconds and remove the catheter slowly from the uterus.
- NOTE:** Please note that after using the "ABORT" feature, the treatment can NOT be resumed nor can the device be restarted.

Steps for the procedure:

- 1. Remove the blister tray from the environmental pouch.
- 2. Carefully inspect the blister trays and Tyvek lids for any damages prior to use. Do NOT attempt to use the device if the sterile barrier is damaged. Do NOT use past expiration date.
- 3. a. Move the indicator ring backwards.
b. Open and remove the balloon cap.
- 4. Insert the battery into the device and place the device on the sterile field.
The LiNA Librata™ is designed with 47 minutes of battery life, that covers one procedure only.
NOTE: Do not press the "ON"-button before all preparations have been made (as the heating process will then begin).
- 5. Conduct a pelvic examination to determine the position of the uterus. Particular attention should be made when sounding, dilating and inserting the device into uterus. Insert a speculum into the vagina and grasp with a tenaculum to bring the uterus to a mid-plane position.
 - Using a uterine sound, measure the length of the uterus from the external os to the fundus. Set the indicator ring on the catheter to the measurement obtained from the sounding. (Please note the markings on the LiNA Librata™ catheter are intentionally 0.5cm shorter than indicated. This allows a 0.5cm gap between the balloon and fundus once the catheter is inserted until the indicator ring reaches the cervix).
- NOTE:** Due to the small diameter of the LiNA Librata™ catheter, dilation of the cervix may not be required. In the event that the cervix needs to be dilated, dilate the cervix to 6mm.
- 6. Press the "ON" button. The device will make a short beep sound and the "HEATING" light will start flashing. The heating process takes approximately 10 min.
NOTE: During the heating process the device should remain on the sterile field and should not be hand-held or inserted into the patient.
- 7. Once the fluid is heated to the correct temperature, the "HEATING" light will stop flashing. The device will now make a short beep sound and the "READY" light will illuminate. The device is now ready and LiNA Librata™ treatment can be initiated.
NOTE: When the device is ready, the LiNA Librata™ treating process needs to be started within 30 minutes. If the treating process is not started within this time frame, the fluid can NOT be reheated and the procedure may therefore be aborted.
NOTE: Please note that the device will make a short beep sound every minute until the treatment is started.

- 8. Carefully insert the catheter into the uterus until the indicator ring touches the cervix. If resistance is felt before the indicator ring is touching the cervix, remove the catheter and take a second sounding of the uterus and adjust the indicator ring accordingly.
NOTE: Lubricant is not required when inserting the catheter.
NOTE: Do not use excessive force when inserting the balloon catheter.
- 9. When the device has been inserted correctly into the uterus, the treating process may begin. Press and hold the "TREAT" button for 2 seconds and the "TREATING" light will then begin to flash.
- 10. The treating process takes approximately 2 minutes. During the treating process fluid will automatically be injected and withdrawn from the balloon a total of 5 times. Each cycle will take approx. 25 seconds.
NOTE: Prior the beginning of the treatment process, the device performs a short 5 second pressure test to indicate that the balloon has not passed outside the uterine cavity. If the device suspects that the balloon is outside the uterine cavity the device will abort the procedure. All lights on the device will flash simultaneously.
A hysteroscopy should then be performed to confirm that a perforation has not taken place. If no perforation is revealed during hysteroscopy, a second LiNA Librata™ device may be used.
- 11. When the treating process is done the "FINALIZED" light will illuminate and the device will make a short beep sound.
- 12. Remove the catheter slowly from the uterus.
NOTE: If unexpected resistance is felt during removal of the catheter pull the "EMERGENCY" strap on the back of the device until the word "STOP" appears. Wait 5 seconds and remove the catheter slowly from the uterus.
- 13. Dispose of the device. The battery and the device may now be safely disposed in accordance with local governing ordinances and recycling plans. If desired the battery can be removed from the device and disposed separately.

DECLARATION OF CONFORMITY

with the Medical Device Directive:
We, LiNA Medical ApS, Formervangen 5, 2600 Glostrup, Denmark under the sole responsibility of the manufacturer declare that the catalogue number(s):
Product name(s) REF number(s)
LiNA Librata™ LIB-1
LiNA Librata™ LIB-3
are in conformity with the European Medical Device Directive 93/42/EEC amended by 2007/47/EC, in effect as of 21st of March 2010 and pertinent national statutory order.

CE marking 2016.

LiNA Librata™

REF: LIB-1, LIB-3

DEVICE ALARMS:

During heat up (before the TREAT button is pushed):

Light pattern	Description	Reason for error	User instruction
	The "FINALIZED" and "READY" lights are flashing during heat up.	Device failure for technical reasons only, no patient safety issue is suspected.	1. Return the device for investigation and dispose of the battery. 2. Replace with a new device.

During treatment (after the TREAT button has been pushed):

Light pattern	Description	Reason for error	User instruction
	The "FINALIZED" and "READY" lights are flashing during treatment.	The system has detected temperature or pressure levels outside the preset limits.	1. Wait 8 seconds to let the "ABORT" feature withdraw all the fluid from the balloon. 2. Perform a hysteroscopy to rule out any possible patient injury. 3. Return the device for investigation and dispose of the battery.
	The "TREATING" and "HEATING" lights are flashing.	The uterine cavity is too big to treat. The patient cannot be treated with the LiNA Librata™ device.	
	All lights are flashing.	The balloon was unable to meet a minimal pressure before initiating treatment. The balloon may be outside the uterine cavity.	

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LINA is a registered trademark of LINA Medical in EU, US and CN