
LiNA SeaStar™

(Reusable)

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

LiNA SeaStar™

REF 4000-R-SL, 4000-R-VL, 4000-R-VS, 4000-R-SS, 4000-R-VV

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

Product description:
The LiNA SeaStar™ is a retractor ring with grooves every 2.5 cm on the convex site. The hinges are controlled by finger screws allowing the two ring halves to remain at any desired angle. The device is provided in a broad range of shapes and sizes which can easily be adjusted to fit various anatomical areas, thereby facilitating different surgical applications. The LiNA SeaStar™, together with LiNA Skin Hook™ is a complete system for retracting delicate soft tissue. The flexible design of the system allows easy adjustment and repositioning of the LiNA SeaStar™ during the surgical procedures and provides optimal exposure of the surgical site.

Intended use:
The LiNA SeaStar™ is intended to retract soft tissue from the operative surface.

Patient population:
All groups dedicated to surgical procedures which require optimal exposure.

- Warnings:**
- Device is provided as a non-sterile. Please follow recommended cleaning and disinfection instructions before the first and each consecutive use. Use cleansers and/or disinfectant with a pH range between 4.5 and 12. Follow the manufacturer's instructions regarding dosage, application time and replacement of the solution.
 - Device is provided as a non-sterile therefore it must be sterilized before the first and each consecutive use. Please follow recommended sterilization instructions.
 - The number of times the device is exposed to autoclaving must not exceed 100. If any errors are detected prior to this amount, dispose of the device in accordance with hospital, administrative and/or local government policy.
 - Cleaning, disinfection and sterilization instructions apply exclusively to the reusable retractor rings. The disposable retractor rings by LiNA Medical may not be sterilized.
 - Do not use past expiration date.
 - Do not use if the device is exposed to non-sterile surfaces before procedure.
 - Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
 - No modification of this equipment is allowed.
 - Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA SeaStar™ should not be used for any other purpose than intended function.

- Precautions:**
- Use of the LiNA SeaStar™ requires adequate training and experience in performing surgical procedures. Intended user profile is a trained surgeon or physician.
 - Some accessories used with the LiNA SeaStar™ (e.g. LiNA Skin Hook™) are single use materials. Single use accessories used with the reusable instruments should be discarded after use.
 - Do not use excessive force during the screw tightening.
 - The instruments are to be handled with care, e.g. do not drop them on the floor. Avoid putting loads on the instruments or bending them. In order to avoid damage to the instruments, do not use metal brushes, metal sponges or the like.
 - The LiNA SeaStar is intended for up to 100 autoclave sterilization cycles. Following the use, the LiNA SeaStar is considered to be contaminated medical waste and may pose a potential infection risk. If the product will not be cleaned and reused, dispose of the device in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.
 - Always have a backup device readily available for immediate use.

Instruction for use:
The following instructions apply exclusively to the reusable retractor rings. The disposable retractor rings by LiNA Medical may not be sterilized.

When using for the first time:
Take the instruments out of the packaging and remove the protective wrapping.

Ultrasonic cleaning:
Open and disassemble the instrument. Place instrument parts in the basket avoiding overlap. Add cleaning agent to the water and set the temperature of the solution according to the instructions of the cleaning agent manufacturer. Cleaning with the ultrasonic cleanser at 35 to 40 kHz takes 5 minutes. Rinse instruments with water and afterwards immediately begin with the cleaning and/or disinfecting procedure.

Cleaning and disinfection:
If possible, always choose mechanical cleaning over manual cleaning. Cleaning solutions and disinfectants are to be replaced if they become contaminated and should be replaced at least once per day.

LiNA SeaStar™

REF 4000-R-SL, 4000-R-VL, 4000-R-VS, 4000-R-SS, 4000-R-VV

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

Manual cleaning:
Use synthetic brushes, water jet pistols or lent-free towels and fresh non-adhesive cleaning agents for manual cleaning. Rinse thoroughly with desalinated water. Place instruments in disinfectant and afterwards use a suitable cleaning agent as described above. Rinse the instruments again with desalinated water and then dry them.

Automatic cleaning and disinfecting device:
Chemically clean at 40 to 60 degrees Celsius after the pre-rinse. Residue from the procedure is to be removed carefully in the 2nd rinse phase. Avoid deformations resulting from neutralization agents. Disinfect the instrument parts thermally at 80 to 95 degrees Celsius. Rinse with desalinated water and then let dry.

Drying:
The instruments can be dried in a drying oven or with pressurized air.

Inspection and check:
Allow the instruments to cool down to room temperature and macroscopically check them for residue and contamination. Instruments that are not free from residue must undergo the entire cleaning procedure again. Check the instruments to see if they function as intended. Deformed, damaged or non-functioning instruments may not be used again and must be discarded.

Packaging:
Use packaging material that is suitable for the instruments and which protects against microbiological contamination during storage.

Sterilization:
Steam sterilization of the instruments is mandatory. Packaging and storage material must be adapted to the instruments and the sterilization procedure.

Steam sterilization:
Place instruments into the baskets making sure that the recommended total weight is not exceeded. It is recommended to store the baskets in towels or fleece and in the appropriate containers. The baskets can also be alternatively placed in paper bags. Autoclave between 3 and 18 minutes at 134 degrees Celsius or between 15 and 30 minutes at 121 degrees Celsius. After autoclaving, dry and let cool to room temperature (23 degrees Celsius, +/- 2 degrees).

Description of use:
Place the LiNA SeaStar™ around the initial incision and insert the hooks in directly opposing quadrants with equal tension. As dissection progresses, move the position of the hook into the deeper layers of the fascia. You may then add new hooks or reposition those already in place to achieve full balanced retraction from any angle. To ease closing maintain constant opposing tension.

Storage:
Sterilized instruments must be stored in dry, dust-free and ventilated rooms without corrosive vapors.

Additional information:
Employees must be trained in the sterilization of instruments. The autoclaves must comply with ISO 17665. Washing and disinfectant systems must be approved in accordance with ISO 15833. The packaging material for the sterile equipment must comply with the ISO 11607 standard. The instructions listed above have been tested by the manufacturer and are suitable for sterilizing SeaStar reusable retractor rings. It is the responsibility of the person in charge of sterilization to ensure that only material, equipment and staff which lead to the desired result are deployed. This includes acceptance and regular routine checks of the process.

Reporting:
Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

Caution:
Federal Law (USA) restricts this device to sale by or on the order of a trained physician.

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LiNA is a registered trademark of LiNA Medical in EU, US and CN