
LiNA OperåScope™

Tubing & Drape Kit

Innovation in Gynecology

LiNA Medical AG - Platz 3 - 6039 Root D4 - Switzerland - Tel: +41 78 617 56 56 - e-mail: info@lina-medical.com - www.linamed.com

LiNA Medical Polska Sp. z o.o. - Rolna 8A - Sady - 62-080 Tarnowo Podgórne - Poland - Tel: +48 61 222 21 21 - e-mail: info.poland@lina-medical.com - www.lina-medical.pl


2025-10 FV0279D

LiNA is a registered trademark of LiNA Medical in EU, US and CN

2025-10 FV0279D

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

LiNA OperâScope™ Tubing & Drape Kit

REF OP-TDK-6-I

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Product description:

The LiNA OperâScope™ Tubing & Drape Kit (Tubing & Drape Kit) is delivered as a sterile, single-use device. The Tubing & Drape Kit is designed to be used specifically with the LiNA OperâScope™ hysteroscope. The device consists of a drape with pouch, an inflow tubing, an outflow tubing and a collection bag.

Intended use:

The LiNA OperâScope™ Tubing & Drape Kit is intended for all diagnostic and/or operative hysteroscopy procedures, to manage the fluids used for distending the uterus.

Patient population:

Adolescent and adult females suited for hysteroscopy.

Contraindications:

This product is not intended for use when hysteroscopy is contraindicated.

Contraindications for hysteroscopy include:

- Acute Pelvic Inflammatory Disease (PID)

Hysteroscopy may be contraindicated by the following conditions, depending on their severity or extent:

- Inability to distend the uterus,
- Cervical stenosis,
- Cervical/vaginal infection,
- Uterine bleeding or menses,
- Known pregnancy,
- Invasive carcinoma of the cervix,
- Recent uterine perforation,
- Medical contraindication or intolerance to anesthesia.

Potential complications:

Potential complications of Continuous Flow Hysteroscopy include: uterine perforation/false passage (resulting in possible injury to bowel, bladder, major blood vessels and ureter), cervical laceration, bleeding, vasovagal episodes, pain, infection, uterine adhesions, air embolism and fluid-related complications (hyponatremia, hypothermia, pulmonary edema, cerebral edema, heart failure).

Warnings:

- The LiNA OperâScope™ Tubing & Drape Kit is provided sterile using ethylene oxide sterilization. Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before the procedure.

- For single use only. Do not reuse, reprocess or re-sterilize the components of LiNA OperâScope™ Tubing & Drape Kit. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device result in risk of device malfunction and/or cross contamination due to residual tissue in the tubing or drape.
- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- No modification of this equipment is allowed.
- Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA OperâScope™ Tubing & Drape Kit should not be used for any other purpose than its intended function.
- Do not contaminate the connection between the saline bag and the tubing of the LiNA OperâScope™ Tubing & Drape Kit. Caution must be exercised when connecting the spike to the saline bag.
- To minimize the risk of air embolism, the hysteroscope and inflow tubing should be primed with the distension medium to eliminate air bubbles before inserting the hysteroscope into the uterine cavity.
- The use of normal saline as a distension medium and limiting the infused volume to less than 1000 ml is recommended.
- The lowest possible fluid distension pressure should be used to obtain an adequate view. Reduce inflow – e.g., by partial or full closure of the LiNA OperâScope™ inflow stopcock or by reducing pressure or height of saline bag. Consider doing this routinely on entering the uterine cavity (start 'low' and build up).
- Careful attention should be paid to fluid inflow and outflow. A fluid deficit (difference between the volume of fluid infused into the uterus and the volume of fluid evacuated from the uterus) should not exceed 2500 mL for isotonic solutions (saline). For elderly patients and others with comorbid conditions including compromised cardiovascular systems, a maximum fluid deficit of 1500 mL is recommended. Electrocardiogram changes should prompt discontinuation of the hysteroscopic intervention.
- Always have a backup device readily available for immediate use. If any malfunction should occur that prevents the correct estimation of fluid outflow, the patient should be carefully monitored with the consideration to terminate the procedure if the fluid deficit is thought to have exceeded the recommended maximum values.

LiNA OperâScope™ Tubing & Drape Kit

REF OP-TDK-6-I

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

- Consult the Instructions for Use of the hysteroscope before performing procedures with the LiNA OperâScope™ Tubing & Drape Kit.

Precautions:

- Intrauterine distension can usually be accomplished with pressures in the range of 35 - 75 mmHg. Unless the systemic blood pressure is excessive, it is seldom necessary to use pressures greater than 75-80 mmHg.
- The LiNA OperâScope™ Tubing & Drape Kit should not be connected to fluid management systems/pumps, as compatibility with such systems has not been verified.
- Use of the LiNA OperâScope™ Tubing & Drape Kit requires adequate training and experience in performing hysteroscopy.
- After use, the device is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging as well as the fluids from the saline bag, the collection bag and the drape pouch in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.

Instruction for use:

Ensure the following other devices/equipment are available for the procedure:

- Infusion bag stand with optional pressure cuff,
- Saline bag (distension medium),
- Gynecological bench (examination bench).

1. Throughout the procedure, use sterile techniques to avoid contaminations.
2. Prior to opening the sterile pouch, hang the saline bag in the stand/pressure cuff.
3. Prepare the LiNA OperâScope™ hysteroscope for use by closing both Luers.
4. Remove the inner pouch from the sterilized blister (**Figure 1**).
5. Remove the drape from the inner pouch (**Figure 2a**).
6. Inspect the drape for labels for left & right hand. Unfold the drape by inserting the left & right hand at the labels and pull them away from each other (**Figure 2b-2d**).
7. Place the drape between the bench and the patient's buttock. Let the pouch hang down towards the floor (**Figure 2e**).
8. Remove the collection bag from the inner pouch. Close the drain valve at the bottom and remove the protection cap from the inlet (**Figure 3a**).
9. Remove the OUTflow tubing from the inner pouch. Connect the Luer to the OUT connector on the LiNA OperâScope™ and connect the other end to the collection bag or put it into the drape pouch (**Figure 3b – 3e**).

10. Remove the INflow tubing from the inner pouch.

⚠ Connect the Luer to the IN connector on the LiNA OperâScope™ (**Figure 3f**).

11. Remove the protection cap from the spike and insert the spike into the saline bag (**Figure 3g**). The slide clamp may be used to temporarily stop the flow from the saline bag, if needed.
12. Carefully inspect the saline bag for material fragments.
13. Apply light pressure to the saline bag by **adjusting the hanging height** of the saline bag or by **using a pressure cuff**.
14. Open the inlet Luer on the LiNA OperâScope™ and flush fluid through the tubing and the LiNA OperâScope™ until all air bubbles have cleared from the tubing. Open the outlet Luer.
15. Insert the LiNA OperâScope™ into the patient and perform the procedure.
16. To distend the uterus, adjust the LiNA OperâScope™ outflow stopcock position and the pressure on the saline bag.
17. Throughout the procedure, monitor the fluid deficit to avoid fluid overload (**Figure 4**).
18. When the procedure is complete, remove the LiNA OperâScope™ from the patient and close the inlet Luer, release the pressure from the saline bag. Disconnect the tubings from the LiNA OperâScope™, saline bag and collection bag.
19. Dispose of the fluids from the saline bag, the collection bag, and the drape pouch. If needed, cut the drain plug on the pouch (**Figure 5**).

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

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