
LiNA OperåScope™

Innovation in Gynecology

LiNA OperåScope™

REF OP-201-6-I, OP-RM-1-UK, OP-RM-1-EU, OP-RM-1-AUS

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

- For single use only. Do not reuse, reprocess or re-sterilize the LiNA OperåScope™. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device may result in risk of device malfunction and/or erroneous pathology specimen collection due to residual tissue in the LiNA OperåScope™.
- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- No modification of this equipment is allowed. Do not attempt to alter the shape of the distal cannula.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before the procedure.
- Low light intensity at the lowest available setting may limit visualization and be unsuitable for surgical procedures. Device operators should check that the light intensity range is suitable for the planned procedure prior to use.
- Do not advance the LiNA OperåScope™ should resistance be experienced while in the patient.
- Do not remove the battery during the procedure.
- For use only with the batteries provided with each LiNA OperåScope™ device – batteries are identified as AA Energizer alkaline batteries in a dedicated battery pack for the LiNA OperåScope™.
- The LiNA OperåScope™ is designed with 60 minutes of battery life.
- If any malfunction should occur during use or battery life expires, stop the procedure immediately, and slowly withdraw the LiNA OperåScope™ and replace with a new LiNA OperåScope™.
- Portable RF communication equipment (including peripherals such as antenna cables and external antennas) should not be used within 30 cm (12 inches) of the LiNA OperåScope™, as this could result in the degradation of the LiNA OperåScope™ performance.
- Use of accessories, transducers and cables other than those specified could result in a change in electromagnetic properties resulting in potential improper operation.
- Do not use High Frequency electrosurgical devices with the LiNA OperåScope™, as high-frequency electromagnetic energy causes image interference on video screens.
- Use of the LiNA OperåScope™ equipment adjacent to other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Avoid submersing the LiNA OperåScope™ handle, battery, HDMI port and/or LCD screen in a fluid bath. Should the LiNA OperåScope™ be submerged, ensure the LiNA OperåScope™ is dried off and remains operational before clinical use. If the LiNA OperåScope™ is nonoperational, replace with a new device.
- To mitigate the risk of perforation, only advance/manipulate the LiNA OperåScope™ while viewing a live camera image and illumination.
- Pay attention when the patient has any implantable electronic device (IED). The perioperative team shall consult with the team managing the implanted device before surgery and as needed. In case the team managing the device is not available, consult with the device manufacturer or the interdisciplinary team in your hospital. When the LiNA OperåScope™ is activated it shall be kept more than ½ inch (1 cm) from the IED.
- Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of the procedures performed. The LiNA OperåScope™ should not be used for any other purpose than intended function.
- Strict fluid intake and output surveillance should be maintained. Intrauterine instillation exceeding 1000 ml should be followed with great care to reduce the possibility of fluid overload.
- The use of normal saline as a distending medium, and limiting the infused volume to less than 1000 ml is recommended to decrease the risk of the above complications.
- Careful attention should be paid to fluid inflow and outflow. A fluid deficit (difference between the volume of fluid infused into the uterus and the volume of fluid evacuated from the uterus) should not exceed 2500 ml for isotonic solutions (saline). For elderly patients and others with comorbid conditions including compromised cardiovascular systems, a maximum fluid deficit of 1500 ml is recommended.

Precautions:

- Before each use or after a change of the LiNA OperåScope™ settings, check to ensure the LiNA OperåScope™ view provides a live image (rather than a stored one) and has the correct image orientation.
- If a continuous flush/irrigation during the LiNA OperåScope™ operation is not maintained, the distal tip operating temperature may exceed 41°C (105.8°F).
- Intrauterine distention can usually be accomplished with pressures in the range of 35 - 75 mmHg. Unless the systemic blood pressure is excessive, it is seldom necessary to use pressures greater than 75-80 mmHg.
- While the LiNA OperåScope™ Recording Module is connected to the LiNA OperåScope™, ensure that the power supply unit is easily accessible if power should need to be disconnected.
- Ensure the LiNA OperåScope™ Recording Module and monitor are outside of the patient's reach.
- Always have a backup device readily available for immediate use. If any malfunction should occur during use, stop the procedure immediately, and slowly withdraw the LiNA OperåScope™ and replace with a new device.
- Use of the LiNA OperåScope™ requires adequate training and experience in performing hysteroscopy and cystoscopy.

LiNA OperåScope™

REF OP-201-6-I, OP-RM-1-UK, OP-RM-1-EU, OP-RM-1-AUS

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

- After use, the device, including the battery, is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.
- In premenopausal women with regular menstrual cycles, the optimal timing for diagnostic hysteroscopy is during the follicular phase of the menstrual cycle after menstruation.
- Do not insert or remove the instruments through LiNA OperåScope™ with excessive force, as this may lead to damage of the LiNA OperåScope™ and instruments.

Instruction for use:

Other devices/equipment which may be used during the procedure with the LiNA OperåScope™:


- The LiNA OperåScope™ Recording Module,
- USB stick,
- HDMI cable,
- Commercially available fluid management tubing with a male luer lock,
- Bag with saline (distention medium),
- Commercially available mechanical instruments for diagnostic/therapeutic use,
- External Monitor with CE mark.

External monitor specification

Parameter	Requirements
Power supply	Monitor rated with maximum Input Voltage 240VAC.
Monitor resolution	As minimum 640x480 pixels or higher with refresh rate 60Hz. Recommended 720p (HD) or 1080p (Full HD).
Video input	Minimum 1x HDMI 1.3 or higher.
Electrical safety	Monitor has to be Medical grade (IEC 60601-1) or provides galvanic separation from mains power supply by compliance with one of the following electrical safety standards: IEC 62368-1, IEC 60950-1 or IEC 60065.

- Using sterile technique, remove the LiNA OperåScope™ from the sterilized blister (**Figure 1**).
- Attach the inflow fluid tubing to the inflow port and the outflow tubing to the outflow port (**Figure 2**).
- The LiNA OperåScope™ can be used in the following configurations:
 - As a standalone device, not connected to any external device (Recording Module or External Monitor).
 - Connected directly to an External Monitor which bears the CE mark via HDMI cable for visualization of the image.
 - Connected to the Recording Module via HDMI cable for recording and storing information related to the procedure.
 - Connected to a Recording Module that is already connected to an External Monitor.

Note: Prior to using the LiNA OperåScope™ Recording Module, ensure to affix the wall plug relevant for the local wall sockets and the unit is plugged in and green power light is illuminated.

If picture or video recording is desired, insert a USB memory stick into the LiNA OperåScope™ Recording Module and make sure the HDMI cable is connected from the LiNA OperåScope™. Both the LiNA OperåScope™ handle and Recording Module have picture and video capture buttons. Press and release the picture button to capture a still image. Press and release video button to start recording; press and release the video button again to stop recording (**Figure 3**).

LiNA OperåScope™

REF OP-201-6-I, OP-RM-1-UK, OP-RM-1-EU, OP-RM-1-AUS

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Note: When the LiNA OperaScope™ Recording Module is busy the LED will be flashing blue.

- Remove the battery from the sterilized blister and insert the battery into the LiNA OperaScope™ handle (**Figure 4**).
- The light source at the distal tip of the cannula will illuminate and the camera image should be visible on the onboard LCD. If illumination or camera image is not functional, discontinue setup and open a new LiNA OperaScope™. The cannula is rotated by a turning knob on the handle (**Figure 5**).
- Flush the fluid through the cannula until all air bubbles have cleared from the tubing (**Figure 6**).
- The LiNA OperaScope™ is now ready for use. With fluid flowing, the cannula is inserted through the cervical canal or urethra while viewing the image, as is standard procedure for introduction of the LiNA OperaScope™ (**Figure 7a**). The LiNA OperaScope™ has a central lumen located within the turning knob of the handle (**Figure 7b**) that provides a conduit for delivery of diagnostic and therapeutic instruments/accessories to the uterine or bladder cavity. Mechanical instruments within the following dimensions are compatible with the LiNA OperaScope™ central lumen:
 - Maximum outer diameter: 1.86 mm (5.5 Fr)
 - Minimum useable length: 310 mm
- After the procedure is completed, remove the LiNA OperaScope™ (**Figure 8**).
- Remove the battery pack and dispose of the battery. The device may now be safely disposed. After use, the device, including the battery, is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances (**Figure 9**).
-  Remove USB stick from the LiNA OperaScope™ Recording Module and transfer the recordings to a computer using a USB port (**Figure 10**).
- After each patient use, power off and unplug the LiNA OperaScope™ Recording Module and disinfect it according to post procedure cleaning section (**Figure 11**).

Technical specification for the LiNA OperaScope™

Nominal input voltage: 3V DC $\pm 10\%$
Nominal Current: 330mA $\pm 20\%$
Battery type: 2pcs AA 1.5V Energizer Alkaline (Non-rechargeable)
Operating time: ≤ 60 minutes
Camera resolution: 400x400 pixel
Diagonal size of LCD display: 4.3 Inch
Minimum usable length: 240mm
Pre-fixed tip angle: $20^\circ \pm 5^\circ$

Shelf Life Conditions

The shelf life for the LiNA OperaScope™ is 36 months. The test units were then subjected to simulated shipment conditions per ASTM D4169, Cycle 13, assurance level 1. All units were inspected and subjected to design verification testing and met all acceptance criteria.

Environmental conditions during transportation and storage

Temperature: 0°C to $+55^\circ\text{C}$
Air humidity: 10% to 90%
Atmospheric pressure: 70 kPa to 106 kPa

Environmental conditions during use:

Temperature: $+10^\circ\text{C}$ to $+40^\circ\text{C}$
Air humidity: 35% to 75%
Atmospheric pressure: 70 kPa to 106 kPa

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

CE marking: 2020.

LiNA OperåScope™

REF OP-201-6-I, OP-RM-1-UK, OP-RM-1-EU, OP-RM-1-AUS

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Guidance and manufactures declaration for electromagnetic immunity				
Emission standard	Type	Compliance level		Specific guidance for application environment
EN/IEC61000-4-2	Electrostatic discharge (ESD)	± 8kV contact Up to ± 15 kV air		The device is suitable for use with mains not connected to the public low-voltage network supplying domestic environment. Portable and mobile radio devices, including their wires, should not be used closer to the unit than 30 cm.
EN/IEC61000-4-4	Electrical fast transients and bursts	± 2 kV power line ± 1 kV signal lines (100 kHz repetition)		
EN/IEC61000-4-5	Surge	± 0.5, ±1kV line to line		
EN/IEC61000-4-11	Dips, interruptions and variations on power supply input	0%, 0.5 cycles @ 0; 45; 90; 135; 180; 225; 270; 315o 0%, 1 cycle 70%, 25/30 cycles Interruptions: 0%, 250/300 cycles		
EN/IEC61000-4-8	Power frequency magnetic field	30 A/m at 50/60 Hz		
EN/IEC 61000-4-6	Conducted disturbances, induced by RF fields	150 kHz to 80 MHz 3 V r.m.s before modulation (6 V ISM bands)		
EN/IEC61000-4-3	Radiated RF EM fields	3 V/m 80 MHz to 2.7 GHz, 80% AM @ 1 kHz		
EN/IEC61000-4-3	Proximity fields from RF wireless communication equipment	Test freq. MHz	V/m	
		385	27	
		450	28	
		710, 745, 780	9	
		810, 870, 030	28	
		1720, 1845, 1970	28	
		2450	28	
		5240, 5500, 5785	9	

Guidance and manufactures declaration for electromagnetic emissions		
Emission standard	Type	Compliance to standard
CISPR 11	RF emissions	Yes, Group 1, Class A
Guidance for application Environment		
The OperåScope™ device is suitable for use in professional healthcare facilities in environments not connected to the public low-voltage network.		
The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential ENVIRONMENT (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.		

LiNA OperåScope™

REF OP-201-6-I, OP-RM-1-UK, OP-RM-1-EU, OP-RM-1-AUS

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

LiNA OperåScope™ Recording Module

REF OP-RM-1-UK, OP-RM-1-EU, OP-RM-1-AUS

Intended Use:

The LiNA OperåScope™ Recording Module is intended for recording of video or pictures to a standard USB drive.

The LiNA OperåScope™ Recording Module power supply is prepared according to the above bullet 3 and 11 in Instructions for use. The different ref. numbers refer to the relevant socket.

- OP-RM-1-EU (for use in most EU countries – socket type C)
- OP-RM-1-UK (for use in UK and other countries using socket type G)
- OP-RM-1-AUS (for use in Australia and other countries using socket type I).

Power supply information

Nominal input voltage:	100-240 V AC ±10%
Nominal input frequency:	50-60 Hz
Nominal input current:	0.16-0.08 Arms @ max load
Nominal output voltage Uout:	9 V DC ±5%
Nominal output current Iout:	800 mA
AC adapter provided:	Friwo Fox 9 V DC/800 mA

The following accessories are included in the LiNA OperåScope™ Recording Module:

Description	Pcs
HDMI-DVI conversion adaptor:	1
HDMI cable 3m:	1
USB memory stick 16GB:	1
AC adaptor:	1
Hook and loop fastener:	2 sets
Wall hook:	1
Zip ties 20mm:	2
HDMI terminal cover:	2

Download of Video/Still Imagery Files Post Procedure

Video and still imagery can be downloaded from the LiNA OperåScope™ via the LiNA OperåScope™ Recording Module USB and imported into the patient's medical record. Since there is no patient identifying information on the video/still imagery file, the downloaded file should be uploaded and saved to the patient's medical record immediately post procedure ensuring the proper link between the file and the patient.

Playback of recorded Video/Still images on computer

Step 1

Insert the USB-stick into the computer USB port.

Step 2

Download and install some supplement codec software for Windows Media Player such as Media Player Codec Pack.

Step 3

Click "Start" in Windows and select the folder in which your MPEG 4 file is located, such as My Videos or My Downloads.

Step 4

Right-click on the video and select "Open With." Choose "Windows Media Player" from the available options to begin playing your video/Images.

LiNA OperåScope™

REF OP-201-6-I, OP-RM-1-UK, OP-RM-1-EU, OP-RM-1-AUS

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Post Procedure Cleaning

The LiNA OperåScope™ Recording Module is intended to be used outside the sterile field and away from the patient. While the LiNA OperåScope™ Recording Module is not intended to come in direct contact with patient, the device can be exposed to body fluids and/or blood as it is used in the operating room during patient treatment.

After each patient use, the following cleaning/disinfecting instructions needs to be completed to ensure proper cleaning of the external components.

1. Power off and unplug the LiNA OperåScope™ Recording Module.
2. Place the LiNA OperåScope™ Recording Module and power cord on a clean, dry surface.
3. Using a cleaning/disinfecting wipe, completely wipe the top, bottom and all 4 sides of the LiNA OperåScope™ Recording Module.
4. After wiping, visually inspect the exterior for residual bodily fluids or dirt and repeat wipe.
5. Let the LiNA OperåScope™ Recording Module sit for at least 3 minutes.
6. Using a 70% Isopropyl swab, thoroughly wipe the seams of the Module ensuring the Module is free of bodily fluid.
7. Once the Module is visually clean, place the device in a storage area for future use.

Precautions related with cleaning:

- Do not submerge the LiNA OperåScope™ Recording Module;
- During cleaning/ disinfecting, orient the LiNA OperåScope™ Recording Module position to minimize liquid from coming in contact with the power, USB, scope and monitor connectors;
- The LiNA OperåScope™ Recording Module case should not be opened for cleaning/ disinfection.

Warnings related with cleaning:

- No modifications of this equipment is allowed.

Safe Disposal

At the end of the product life, clean and disinfect the LiNA OperåScope™ Recording Module as per the instructions in the LiNA OperåScope™ Recording Module Post Procedure Cleaning section of the IFU.

If the USB flash drive is to be disposed of as well, it is recommended to delete the content of the drive.

Then, dispose of the LiNA OperåScope™ Recording Module in accordance with applicable laws, administrative, and/or local government policy and hospital procedures.

General information:

OperåScope™ Recording Module battery durability is 7 years.

Environmental conditions during transportation and storage

Temperature:	0°C to +55°C
Air humidity:	10% to 90%
Atmospheric pressure:	70 kPa to 106 kPa

Environmental conditions during use:

Temperature:	+10°C to +40°C
Air humidity:	35% to 75%
Atmospheric pressure:	70 kPa to 106 kPa

The LiNA OperåScope™ complies with standards for EMC: Medical Electrical Equipment - PART 1-2: General requirements for basic safety and essential performance.

- Note: the emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals/clinics (CISPR 11 Class A). The equipment is not indicated for use outside a hospital or clinic setting.
- Immunity according to Professional healthcare facility environment.

LiNA OperåScope™

REF OP-201-6-I, OP-RM-1-UK, OP-RM-1-EU, OP-RM-1-AUS

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Cybersecurity

The LiNA OperaScope™ Recording Module includes a USB port for storing and displaying timestamped images and videos of the current patient. No personally identifiable data is stored.

Use only the USB drive provided by LiNA Medical for image and video export. It is recommended to scan the USB drive regularly using antimalware or antivirus software.

The LiNA OperaScope™ and the LiNA OperaScope™ Recording Module are equipped with an HDMI port for video output only, intended for connection to a video display.

Software updates for the LiNA OperaScope™ and the LiNA OperaScope™ Recording Module can only be performed by returning the device to LiNA Medical. Do not attempt to make any modifications to the devices.

The Manufacturer Disclosure Statement for Medical Device Security (MDS2) and the Software Bill of Materials (SBOM) are available upon request.

Cybersecurity support will terminate concurrently with the device end of life. Should you require additional information concerning the end of life of the product, please contact your local distributor.

Continued use of the device beyond its end-of-life may increase cybersecurity risks and is not recommended.

Configuration of compatible firmware

LiNA OperaScope™ FW version	LiNA OperaScope™ Recording Module FW version
V5.0	V3.0, V4.0
V4.1	V3.0

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LiNA is a registered trademark of LiNA Medical in EU, US and CN