

LiNA OperåScope™

Scissors

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

LiNA OperâScope™ Scissors

REF OP-SCI-3-I

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Product description:

The LiNA OperâScope™ Scissors (Scissors) are delivered as a sterile, single-use device. The Scissors are designed to be used with the LiNA OperâScope™ hysteroscope. The device consists of handles, a shaft and the blade system (scissors) at the distal end. The Scissors are compatible with hysteroscopes with a straight working channel of min Ø1.87 mm diameter and a max length of 310mm.

Intended use:

The LiNA OperâScope™ Scissors are intended for use in hysteroscopy procedures to cut soft tissue such as endometrial polyps, intrauterine adhesions or septa.

Patient population:

Adolescent and adult females suited for hysteroscopy.

Contraindications:

This product is not intended for use when hysteroscopy is contraindicated.

Contraindications for hysteroscopy include:

- Acute Pelvic Inflammatory Disease (PID)

Hysteroscopy may be contraindicated by the following conditions, depending on their severity or extent:

- Inability to distend the uterus,
- Cervical stenosis,
- Cervical/vaginal infection,
- Uterine bleeding or menses,
- Known pregnancy,
- Invasive carcinoma of the cervix,
- Recent uterine perforation,
- Medical contraindication or intolerance to anesthesia.

Potential complications:

Potential complications of Continuous Flow Hysteroscopy include: uterine perforation/false passage (resulting in possible injury to bowel, bladder, major blood vessels and ureter), cervical laceration, hemorrhage, vasovagal episodes, pain, infection, uterine adhesions, air embolism and fluid-related complications (hyponatremia, hypothermia, pulmonary edema, cerebral edema, heart failure).

Warnings:

- The LiNA OperâScope™ Scissors are provided sterile using ethylene oxide sterilization. Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before procedure.
- For single use only. Do not reuse, reprocess or re-sterilize the LiNA OperâScope™ Scissors. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device result in risk of device malfunction and/or erroneous pathology specimen collection due to residual tissue in the LiNA OperâScope™ Scissors.
- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- No modification of this equipment is allowed.
- Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA OperâScope™ Scissors should not be used for any other purpose than their intended function.
- To mitigate the risk of perforation, only advance or manipulate the LiNA OperâScope™ Scissors while viewing a LiNA OperâScope™ live camera image, allowing observation of the uterine cavity.
- Do not advance the device should resistance be experienced while in the patient.
- Consult the Instructions for Use of the hysteroscope before performing procedures with the LiNA OperâScope™ Scissors.

Precautions:

- Use of the LiNA OperâScope™ Scissors requires adequate training and experience in performing hysteroscopy.
- Always have a backup device readily available for immediate use. If any

LiNA OperâScope™ Scissors

REF OP-SCI-3-I

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

malfunction should occur during use, stop the procedure immediately, and slowly withdraw the LiNA OperâScope™ Scissors and replace with a new device.

- After use, the device is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.
- The LiNA OperâScope™ Scissors are intended only as an adjunct in assessing patient condition. It must be used in conjunction with clinical signs and symptoms.

Instruction for use:

1. Inspect the label to ensure that the correct device is chosen (**Figure 1**).
2. Using sterile technique, remove the device from the sterilized blister (**Figure 2**).
3. Remove the protection cap (**Figure 3a**).
4. Inspect the Scissors blades for any obvious damage.
5. Check the open/close function of the Scissors by actuating the handles. Opening and closing the blade system is done by pushing or pulling the Handles from/ towards each other (**Figure 3b & 3c**).
6. Close the Scissors (**Figure 4a**) before inserting the device into the LiNA OperâScope™ working channel (**Figure 4b**).
7. Slowly insert the Scissors into the LiNA OperâScope™ working channel until the blades are visualized on the LiNA OperâScope™ live camera image. Some resistance will be felt when the blades system passes through the LiNA OperâScope™ pre-curved tip. Take care not to bend the Scissors shaft.
8. Maneuver the Scissors towards the targeted site. Open the Scissors by activating the handles (**Figure 4c**). Advance the opened Scissors over to the targeted site, rotate the Scissors, if needed, (**Figure 4d**) and close the Scissors. Use only enough force to cut the tissue. Over-exertion could cause the Scissors to fail.

9. Continue to apply gentle force on the handles (**Figure 4e**) as the Scissors are extracted from the working channel (**Figure 4f**).

10. Steps 4– 9 may be repeated if multiple cuts are needed within a single procedure.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

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