

LiNA

Versascope™ Sheath

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

LiNA Versascope™ Sheath

REF GMS805

Device description:

The LiNA Versascope™ Sheath GMS805 is a curved access sheath of small diameter intended to be used only with the compatible hysteroscope AlphaScope Hysteroscope (GMS40A). The LiNA Versascope™ Sheath is a single use, sterile device comprised of both metal and plastic components. The device incorporates a channel for insufflation accessed through the attached tubing and connector. Also incorporated is an expandable channel for insertion of hysteroscopic instruments which provides continuous flow where fluids media is utilized for distention.

Intended use:

The LiNA Versascope™ Sheath is used to establish and maintain distention in the uterus and provide access to the uterine cavity for the compatible hysteroscope and other hysteroscopic instruments during diagnostic and operative hysteroscopic procedures.

Patient population:

Adult females undergoing diagnostic or therapeutic hysteroscopy with the compatible hysteroscope.

Contraindications:

This product is not intended for use when hysteroscopic procedures are contraindicated.

Potential complications:

When using a fluid distention medium, strict fluid intake and output monitoring is required. Excessive intravasation of distention fluid can lead to fluid overload. Potential fluid-related complications of Continuous Flow Hysteroscopy are:

- Hyponatremia,
- Hypothermia,
- Pulmonary edema,
- Heart Failure,
- Electrolyte Imbalances,
- Cerebral edema,
- Encephalopathy.

Other potential complications of Hysteroscopy include:

- Injury of the uterine wall/uterine perforation resulting in possible injury to bowel, bladder, major blood vessels and ureter,
- Cervix injury,

- Bleeding,
- Gas embolism,
- Pain,
- Intrauterine Adhesions,
- Infections,
- Vasovagal reactions.

Warnings:

- The LiNA Versascope™ Sheath is supplied sterile using an ethylene oxide. Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged.
- Hysteroscopy procedures performed during pregnancy can result in miscarriage. Suspicion of pregnancy should suggest a pregnancy test prior to the performance of diagnostic hysteroscopy.
- When hysteroscopic instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- Only fluid distension media can be used with the Sheath. LiNA Versascope™ Sheath is not applicable for hysteroscopy using CO₂ as distension medium.
- Strict fluid intake and output monitoring is required. Excessive intravasation of distension fluid can lead to fluid overload.
- Failure to follow all applicable instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA Versascope™ Sheath should not be used for any other purpose than intended function.
- Refer to appropriate electro-surgical system user manual indications and instructions to ensure that all safety precautions are taken.
- A thorough understanding of the principles and techniques involved in laser, electro-surgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device and other medical instruments. Ensure that insulation or grounding is not compromised. Do not immerse electro-surgical instruments in liquid, unless the instruments are specifically designed and labeled to function in liquid.

en

LiNA Versascope™ Sheath

REF GMS805

- For single use only. Do not reuse, reprocess or re-sterilize the LiNA Versascope™ Sheath. Any reprocessing may impede the functions of the device. Reusing single use devices may also increase the risk of cross contamination which may lead to infection or transmission of bloodborne pathogens to patients and users. Attempts to clean the device results in risk of device malfunction.
- Do not use past expiration date.
- Do not use device if it has been exposed to non-sterile surfaces before procedures.
- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- No modification of this equipment is allowed.

Precautions:

- After use, the device is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.
- Only persons having adequate training and familiarity with hysteroscopy should perform hysteroscopic procedures. Consult the IFU of the compatible hysteroscope as well as medical literature relative to techniques, complications, and hazards prior to performance of any hysteroscopic procedure.
- Vaginal ultrasonography before hysteroscopy may identify clinical conditions that will alter patient management. Intrauterine distension can usually be accomplished with pressures in the range of 35–75 mm Hg.
- Unless the systemic blood pressure is excessive, it is seldom necessary to use pressures greater than 75–80 mm Hg.
- A fluid deficit (difference between the volume of fluid infused into the uterus and the volume of fluid evacuated from the uterus) should not exceed 1000 mL for hypotonic (glycine, mannitol, sorbitol) and 2500 mL for isotonic solutions (Umrinikar S, Clark TJ, Saridogan E, et al. British Society for Gynaecological Endoscopy /European Society for Gynaecological Endoscopy Guideline Development Group for Management of Fluid Distension Media in

Operative Hysteroscopy. BSGE/ESGE guideline on management of fluid distension media in operative hysteroscopy. Gynecol Surg. 2016; 13(4): 289–303),

- In case of elderly women with renal and cardiovascular problems, the upper fluid deficit levels should be 750 mL and 1500 mL, respectively. Therefore, careful monitoring of the fluid deficit during hysteroscopy is necessary (AAGL Practice Report: Practice Guidelines for the Management of Hysteroscopic Distending Media. J Minim Invasive Gynecol. 2013; 20(2): 137–148).


How supplied:

The LiNA Versascope™ Sheath is supplied sterile in a Tyvek pouch for single use only. The device will remain sterile until the expiration date as long as the package is unopened and undamaged.

Device can be used in office, hospitals and ambulatory surgical suite procedure room by trained professionals.

Instructions for use:

1. Visually inspect for breaches of sterile barrier system prior to use. Open the package and remove the LiNA Versascope™ Sheath using aseptic handling and transfer methods. Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.

2. Remove the protective cover from the introducer tip.
3. Attach an irrigation line to the insufflation tube of the introducer and an outflow line to the outflow line connector.
4. Place the sterile AlphaScope Hysteroscope through the scope port of the Sheath until the Scope snaps in place.

LiNA Versascope™ Sheath

REF GMS805

5. Attach the hysteroscope to the light source and camera, if desired.
6. With the clamp on the insufflation tube open, and outflow tube clamp closed, begin flow of distension media prior to introducing the LiNA Versascope™ Sheath System into the cervical canal. Slow introduction into the cervical canal will allow the distension media to enhance ease of entry.
7. With the clamp on the insufflation tube open, and the outflow tube clamp closed, insufflate to desired level of pressure and volume or until adequate field of view is obtained. If outflow is desired, open the roller clamp on the outflow tube and insert either a Hysteroscopic Instrument or a LiNA Versascope™ Outflow Cannula. Regulate desired distension and outflow by adjusting the roller clamp on the outflow tube. **WARNING:** Failure to open clamp on the outflow tube and/or failure to use the outflow cannula, may result in overdistension of the uterus or excessive intravasation of fluid.
8. Before starting therapeutic procedure make sure instruments are in closed position when inserting.
9. At end of procedure, ensure any instruments within the instrument port are in the closed position or have been withdrawn from the instrument port prior to withdrawing the Hysteroscopic Instrument System. Instruments should be reprocessed and/or disposed of according to instructions in the instrument IFU. Close inflow clamp and open roller clamp to remove distension media from uterine. **WARNING:** Failure to close the instruments may result in damage to the sheath and potentially pieces being left behind in patient.
10. AlphaScope should be reprocessed and disposed according to instructions in the hysteroscope IFU. LiNA Versascope™ Sheath and its packaging should be disposed after use in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.

Compatibility with other devices:

- Compatibility between hysteroscopic instrument and accessories from different manufacturers, when employed together, shall be verified prior to initiating the procedure.
- There is no guarantee that instruments selected solely using maximum insertion portion width and working length will be compatible in combination.
- There is no guarantee that instruments selected solely using this minimum instrument channel width will be compatible in the combination.
- LiNA Versascope™ Sheath can be used with 6 Fr instruments.

Storage:

No special storage conditions required. Do not use after expiration date.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

CE marking: 2022.

Country of origin: Poland.

*AlphaScope is a registered trademark of GIMMI® GmbH

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LINA is a registered trademark of LINA Medical in EU, US and CN
*AlphaScope é uma marca comercial registada da GIMMI® GmbH