


LiNA OperåScope™ Lasso

For hysteroscopic polypectomies



Innovation in Gynecology

LiNA OperåScope™ Lasso

For hysteroscopic polypectomies

The LiNA OperåScope™ Lasso is a single use mechanical instrument, specifically designed for removing intrauterine polyps.

Versatile

- Available in two sizes


Ready to use, out of the box

- Sterile, single-use

No energy

- Mechanical instrument

Compatible with hysteroscopes with a straight working channel of minimum Ø1.87 mm diameter and a maximum length of 310 mm.



Ordering Information:

LiNA OperåScope™ and LiNA OperåScope™ Lasso **CE** 0459

| Product code | Description | UOM |
|--------------|--|------------|
| OP-201-6-I | LiNA OperåScope™ with HDMI cable & LCD (sterile, single-use) | Box of 6ea |
| OP-LAS10-3-I | LiNA OperåScope™ Lasso 10mm Basket (sterile, single-use) | Box of 3ea |
| OP-LAS16-3-I | LiNA OperåScope™ Lasso 16mm Basket (sterile, single-use) | Box of 3ea |

Ordering Information:

LiNA Versascope™ **CE** 2274

| Product code | Description | UOM |
|--------------|-------------------------|------------|
| GMS805 | LiNA Versascope™ Sheath | Box of 5ea |

Distributor

©2024 LiNA Medical ApS. All rights reserved. All products may not be available for sale in all geographies. Please contact your local LiNA Medical representative for information about products available in your area. This information is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider. Before using any LiNA product, always refer to the IFU included in our packaging for complete indications for use, contraindications, warnings, precautions, and instructions.

LiNA is a registered trademark of LiNA Medical in EU, US and CN.

LiNA Medical ApS

• Formervangen 5
DK-2600 Glostrup
Denmark

• Tel.: +45 43 29 66 66
Fax: +45 43 29 66 99
info@lina-medical.com

• www.lina-medical.com
CVR: 10054974

EN ISO 13485:2016
Quality Management
System

