

Cystoscopy without
capital cost or
complexity



LiNA OperåScope™

Single-use portable cystoscope



LiNA OperåScope™

Single-use portable cystoscope

The LiNA OperåScope™ single-use portable cystoscope is the first and only fully disposable operative cystoscope. Designed for use in visualization of the female urethra and the bladder cavity during diagnostic and therapeutic cystoscopic procedures.

- **Convenient solution** without capital cost or complexity
- **Single-use** design is ready for use out of the box
- **360° steerable, deflected tip** allows for precise control and maneuverability
- **Slim 4.2mm catheter** designed for improved patient tolerability
- **5.5 french operative channel** for therapeutic use

Ordering information:

Product Code	Description	UOM
OP-201-6-I	LiNA OperåScope™ with HDMI cable & LCD	Box of 6
OP-RM-1-UK	LiNA OperåScope™ Recording Module – Socket Type G	Each
OP-RM-1-EU	LiNA OperåScope™ Recording Module – Socket Type C	Each
OP-RM-1-AUS	LiNA OperåScope™ Recording Module – Socket Type I	Each

LiNA OperåScope™

CE
2274



LiNA Medical AG - Platz 3 - 6039 Root D4 - Switzerland
Phone: +41 78 617 56 56, e-mail: info@lina-medical.com, www.linamed.com



LiNA Medical Polska Sp. z o.o. - Rolna 8A - Sady - 62-080 Tarnowo Podgórne - Poland
Phone: +48 61 222 21 21, e-mail: info.poland@lina-medical.com, www.lina-medical.pl

Please refer to the Instruction for Use for complete instructions, contraindications, warnings and precautions.

EU MDR 2017/745 device class for REF OP-201-6-I: IIa.

Name and number of the Notified Body: TÜV Nord Polska Sp. z o.o. # 2274.

REF: OP-RM-1-EU, OP-RM-1-UK, OP-RM-1-AUS: class I.

EN ISO 13485:2016

Quality Management System

©2025 LiNA Medical AG. All rights reserved.

LiNA is a registered trademark of LiNA Medical in EU, US and CN.

LiNA OperaScope Cystoscopy Brochure / 2025-10 BV0249B en

